

STERILIZATION CONTAINERS & ACCESSORIES



STERILIZATION UNDER CONTROL

Made of anodized aluminium, features an additional cover over the perforation field. This makes it particularly suitable for longer-storage or external transports. For storage over extended periods, the protective lid reduces the dust exposure of the filters.

sharpline.de

SHARPLINE.




All rights reserved.
No part of this catalog may be reproduced in
any form without the written permission from SHARPLINE



SHARPLINE.

PRODUCT CODE INFORMATION

01 . 12 . 11 G

MODEL NAME	MODEL SIZE	MODEL HEIGHT	LID or LID COLOR
01 Standard Bottom Non Perforated 	12 Half Size 	11 Height 	G Green 
02 Standard Bottom Perforated 	34 Middle Size 	13 Height 	Y Yellow 
03 Plasma Bottom Perforated 	11 Full Size 	15 Height 	R Red 
04 Mini Size Bottom Non Perforated 	35 Extra Long Size 	20 Height 	B Blue 
05 Mini Size Bottom Perforated 	36 Extra Large Size 	26 Height 	K Black 



SHARPLINE

We are interested in connecting with individuals that are trying to innovate our world into a better place and willing to share creativity. If we can help in anyway please feel free to email us at info@sharpine.de





We Believe

Developing a creative mindset leads to innovative technologies that can change the world. Our passion is to bring products to market that can enhance lives. One of the key ingredients to our success is having 25 years of producing and marketing experience.



SUITABILITY OF OUR STERILE CONTAINERS FOR STEAM STERILIZATION PROCEDURES ¹⁾
as per DIN 58953 part 9

Sterilization container	Packing		Steam sterilization procedure acc. to DIN 58946, part 1				
	Description	Symbolic representation and packing system	Gravitation procedure	Pre-vacuum procedure	Fractional vacuum procedure	Fractional flow procedure	Steam injection procedure
with filters (StB-F)	Textiles in sterilization container with filter in lid and bottom acc. to DIN 58952 part 1		+	+	+	+	+
			○	○	+	+	+
	Textiles in sterilization container with filter in lid, bottom not perforated acc. to DIN 58952 part 1		-	-	+	+	+
			-	-	+	+	○
	Textiles in sterilization container with filter in bottom lid not perforated acc. to DIN 58952 part 1		○	+	+	+	+
			-	○	+	+	+
	Instruments in wire basket acc. to DIN 58952 part III, in sterilization container with filter in lid and bottom acc. to DIN 58952 part 1		+ ³⁾	+	+	+	+
	Instruments in wire basket acc. to DIN 58952 part III, in sterilization container with filter in lid and bottom not perforated acc. to DIN 58952 part 1		-	+ ²⁾	+ ²⁾	+ ²⁾	+ ²⁾
	Instruments in wire basket acc. to DIN 58952 part III, in sterilization container with filter in bottom acc. to DIN 58952 part 1		○	+	+	+	+
with Valves (StB-V)	Textiles in sterilization container with valves in lid and bottom acc. to DIN 58952 part 1		-	○	+	+	○
			-	-	+	+	○
	Textiles in sterilization container with valves in lid no valves on the bottom		-	-	+	+	-
			-	-	+	+	-
	Instruments in wire basket in sterilization container with valve in lid and bottom		-	+ ²⁾	+ ²⁾	+ ²⁾	+ ²⁾
	Instruments in wire basket in sterilization container with valve in lid and not perforated		-	+ ²⁾	+ ²⁾	+ ²⁾	+ ²⁾

Explanation of symbols:

⊕ applicable

⊖ not applicable

○ limited application

- 1) In order to avoid mistakes with regards to service, operation and/or when loading the sterilizer, please refer to table 1 for correct or wrong combination of sterilization containers for standardized steam sterilization procedures.
- 2) Only sterilization containers guaranteeing evaporation or discharge of condensate are suitable above.
- 3) Gravitation procedures having no drying device do not guarantee a sufficient after-drying effect.

INDEX

Standard Model Half Size	8 - 9
Standard Model Middle Size	10 - 11
Standard Model Full Size	12 - 13
Standard Model Extra Long Size	14 - 15
Standard Model Extra Large Size	16 - 17
Plasma Model (Half, Middle, Full Size)	18 - 19
Small Model Mini Size	20 - 21
Small Model Quarter Size (Dental Model)	22 - 23 - 24 - 25
Small Model Implant Size	26 - 27
Small Model Endoscopic Size	28 - 29 - 30 - 31
Small Model Scopy Size	32 - 33
Wire Baskets Half Size	34 - 35
Wire Baskets Middle Size	36
Wire Baskets Full Size	37
Perforated Trays Half Size	38 - 39
Perforated Trays Middle Size	40
Perforated Trays Full Size	41
Perforated Trays Extra Long Size	42
Perforated Trays Extra Large Size	43
Accessories	44 - 51

Standard Model **Half Size**


Superior construction of anodized aluminium
that strenghtens 15%-20% more
than containers of other brands.



SHARPLINE.

Half Size

STERILIZATION CONTAINERS



Bottom
Non Perforated



Outside Dimension
310 x 280 x 110 mm



Outside Dimension
310 x 280 x 130 mm




Outside Dimension
310 x 280 x 150 mm



Outside Dimension
310 x 280 x 200 mm

Safety Lids



Bottom
Perforated



Outside Dimension
310 x 280 x 110 mm



Outside Dimension
310 x 280 x 130 mm



Outside Dimension
310 x 280 x 150 mm



Outside Dimension
310 x 280 x 200 mm



Outside Dimension
310 x 280 x 250 mm



Standard Model Middle Size


SHARPLINE standard model container systems are designed to perfect steam flow and penetration for loads in accordance to the standards and can be used with single use paper filter or reusable PTFE (Teflon) filters.



SHARPLINE.

Middle Size

STERILIZATION CONTAINERS




Bottom
Non Perforated



Outside Dimension
475 x 285 x 110 mm





Outside Dimension
475 x 285 x 130 mm



Outside Dimension
475 x 285 x 155 mm

Safety Lids




Bottom
Perforated



Outside Dimension
475 x 285 x 110 mm



Outside Dimension
475 x 285 x 130 mm



Outside Dimension
475 x 285 x 155 mm



Standard Model Full Size


Warrants sterilization of enclosed contents
and maintains sterilization during
transportation and storage until opened.



SHARPLINE.

Full Size

STERILIZATION CONTAINERS




Bottom Non Perforated



Outside Dimension
600 x 285 x 110 mm



Outside Dimension
600 x 285 x 135 mm




Outside Dimension
600 x 285 x 155 mm



Outside Dimension
600 x 285 x 205 mm

Safety Lids



Bottom Perforated




Outside Dimension
600 x 285 x 110 mm




Outside Dimension
600 x 285 x 135 mm



Outside Dimension
600 x 285 x 155 mm



Outside Dimension
600 x 285 x 205 mm



Outside Dimension
600 x 285 x 255 mm



Standard Model Extra Long Size

Mostly for endoscopic and long (lengthy) instruments.


Caution : Extra Long Sterilization Container can be sterilized only in sterilizers of an appropriate chamber size (at least 726 mm effective depth).



SHARPLINE.

Extra Long Size

STERILIZATION CONTAINERS




Bottom
Non Perforated




Outside Dimension
725 x 285 x 140 mm




Outside Dimension
725 x 285 x 205 mm



Bottom
Perforated



Outside Dimension
725 x 285 x 140 mm



Outside Dimension
725 x 285 x 205 mm

Safety Lids




01.35
01.35.Y
01.35.R
01.35.B
01.35.G
01.35.K



Standard Model Extra Large Size

For Bookwalter Retractor, Craniotomy Set,
Harrington Rod Instrument / Implants, Omni
Retractor, Upper Hand Retractor.




WIDE BODY

sterilization containers


SHARPLINE.

Extra Large Size

STERILIZATION CONTAINERS




Bottom
Non Perforated




Outside Dimension
600 x 372 x 120 mm



Outside Dimension
600 x 372 x 185 mm



Bottom
Perforated



Outside Dimension
600 x 372 x 120 mm



Outside Dimension
600 x 372 x 185 mm

Safety Lids



Plasma Model Low Temperature

Low Temperature
sterilization containers

Also For STERRAD®


STERRAD is a registered trademark owned
by Johnson & Johnson Company



SHARPLINE.

Low Temperature

ALL SIZES
STERILIZATION CONTAINERS



**Bottom
Perforated**



For other models : Wide-Body,
Quarter Size, Mini and Extra Long
Containers please contact us.

Half Size

Code 03.12.11

Outside Dimension
310 x 280 x 110 mm

Code 03.12.13

Outside Dimension
310 x 280 x 130 mm

Code 03.12.15

Outside Dimension
310 x 280 x 150 mm

Middle size

Code 03.34.11

Outside Dimension
475 x 285 x 110 mm

Code 03.34.13

Outside Dimension
475 x 285 x 130 mm

Code 03.34.15

Outside Dimension
475 x 25 x 155 mm

Full size

Code 03.11.11

Outside Dimension
600 x 285 x 110 mm

Code 03.11.13

Outside Dimension
600 x 285 x 135 mm

Code 03.11.15

Outside Dimension
600 x 285 x 155 mm



Small Models **Mini Size**

Mostly Used For:

Meatotomy Set, Newborn Circumcision Set,
Tuboplasty Set, Nasal Local Try,
Eye Surgical Instruments.




SHARPLINE.

Mini Size

STERILIZATION CONTAINERS

Bottom

Non Perforated



	04.14.04
	04.14.04Y
	04.14.04R
	04.14.04B
	04.14.04G
	04.14.04K

Outside Dimension
315 x 142 x 45 mm

	04.14.07
	04.14.07Y
	04.14.07R
	04.14.07B
	04.14.07G
	04.14.07K


Outside Dimension
315 x 142 x 70 mm

	04.14.10
	04.14.10Y
	04.14.10R
	04.14.10B
	04.14.10G
	04.14.10K

Outside Dimension
315 x 142 x 100 mm

Bottom

Perforated



	05.14.04
	05.14.04Y
	05.14.04R
	05.14.04B
	05.14.04G
	05.14.04K

Outside Dimension
315 x 142 x 45 mm

	05.14.07
	05.14.07Y
	05.14.07R
	05.14.07B
	05.14.07G
	05.14.07K

Outside Dimension
315 x 142 x 70 mm

	05.14.10
	05.14.10Y
	05.14.10R
	05.14.10B
	05.14.10G
	05.14.10K

Outside Dimension
315 x 142 x 100 mm



Small Models Quarter/Dental Size

Mostly Used For:

Ear Picks, I & T Tray, Circumcism Set,
Myringotomy Instrument Set.




SHARPLINE.

Quarter Size

STERILIZATION CONTAINERS

Bottom

Non Perforated



Outside Dimension
325 x 192 x 55 mm




Outside Dimension
325 x 190 x 70 mm



Outside Dimension
325 x 190 x 135 mm

Bottom

Perforated



Outside Dimension
325 x 192 x 55 mm



Outside Dimension
325 x 190 x 70 mm



Outside Dimension
325 x 190 x 135 mm



Small Models

1/2 Quarter/Dental Size

Mostly Used For:

Dental small sets, assorted small surgical instruments.





SHARPLINE.

1/2 Quarter Size

STERILIZATION CONTAINERS

Bottom

Non Perforated



06.15.06
06.15.06Y
06.15.06R
06.15.06B
06.15.06G
06.15.06K

Outside Dimension
232 x 174 x 63 mm



06.15.10
06.15.10Y
06.15.10R
06.15.10B
06.15.10G
06.15.10K

Outside Dimension
232 x 174 x 100 mm




06.15.13
06.15.13Y
06.15.13R
06.15.13B
06.15.13G
06.15.13K

Outside Dimension
232 x 174 x 130 mm

Bottom

Perforated



07.15.06
07.15.06Y
07.15.06R
07.15.06B
07.15.06G
07.15.06K

Outside Dimension
232 x 174 x 63 mm



07.15.10
07.15.10Y
07.15.10R
07.15.10B
07.15.10G
07.15.10K

Outside Dimension
232 x 174 x 100 mm



07.15.13
07.15.13Y
07.15.13R
07.15.13B
07.15.13G
07.15.13K

Outside Dimension
232 x 174 x 130 mm



Small Models Implant Size

No tearing anymore. Implant model insures; sealed sterilization, safe transportation and zero contamination.



SHARPLINE.

Implant Size

STERILIZATION CONTAINERS

Bottom

Non Perforated



08.15.07
08.15.07Y
08.15.07R
08.15.07B
08.15.07G
08.15.07K

Outside Dimension
300 x 150 x 75 mm



Bottom

Perforated



09.15.07
09.15.07Y
09.15.07R
09.15.07B
09.15.07G
09.15.07K

Outside Dimension
300 x 150 x 75 mm

Bottom

Non Perforated



08.17.07
08.17.07Y
08.17.07R
08.17.07B
08.17.07G
08.17.07K

Outside Dimension
500 x 170 x 75 mm



Bottom

Perforated



09.17.07
09.17.07Y
09.17.07R
09.17.07B
09.17.07G
09.17.07K

Outside Dimension
500 x 170 x 75 mm

Small Models Endoscopic Size

Endoscopy model provide enough space for all types of endoscopic systems and it comes with perforated lids, non-perforated box and perforated boxes which allowing this system to be used in Plasma Sterilization method.



SHARPLINE

Endoscopic Size

STERILIZATION CONTAINERS



Bottom
Non Perforated



Outside Dimension
535 x 160 x 75 mm



Outside Dimension
535 x 160 x 100 mm



Bottom
Perforated



Outside Dimension
535 x 160 x 75 mm




Outside Dimension
535 x 160 x 100 mm



Basket

Code
21.15.05L
21.15.06L

Outside Dimension
492 x 134 x 55 mm
492 x 134 x 65 mm



Small Models **Plasma Endoscopic**

Endoscopy model provide enough space for all types of endoscopic systems and it comes with perforated lids, non-perforated box and perforated boxes which allowing this system to be used in Plasma Sterilization method.



SHARPLINE.

Plasma Endoscopic

STERILIZATION CONTAINERS



Bottom
Perforated



12.15.07
12.15.07Y
12.15.07R



12.15.07B
12.15.07G
12.15.07K

Outside Dimension: 535 x 160 x 75 mm



12.15.10
12.15.10Y
12.15.10R



12.15.10B
12.15.10G
12.15.10K


Outside Dimension: 535 x 160 x 100 mm



Basket

Code
21.15.05L
21.15.06L

Outside Dimension
492 x 134 x 55 mm
492 x 134 x 65 mm




Small Models Scopy Size

Rigid, Reusable Sterilization container

with silicone holders to maintain

fixation of scopy system during

sterilization and transportaion.




SHARPLINE.

Scopy Size

STERILIZATION CONTAINERS

Bottom

Perforated



13.15.06
13.15.06Y
13.15.06R
13.15.06B
13.15.06G
13.15.06K

Outside Dimension
275 x 70 x 65 mm




13.15.07
1315.07Y
13.15.07R
13.15.07B
13.15.07G
13.15.07K

Outside Dimension
470 x 82 x 74 mm

Bottom

Non Perforated




14.15.06
14.15.06Y
14.15.06R
14.15.06B
14.15.06G
14.15.06K

Outside Dimension
275 x 70 x 65 mm




14.15.07
14.15.07Y
14.15.07R
14.15.07B
14.15.07G
14.15.07K

Outside Dimension
470 x 82 x 74 mm



Standard Models Wire Baskets

SHARPLINE BASKETS for all your process of cleaning the surgical instruments organizing, storing and protecting them even during transportation and sterilization.



SHARPLINE.

Half Size

WIRE BASKETS

Baskets and Perforated Trays are DIN-sized,
so they are compatible with
European-made closed container systems.



Wire Basket

Code	Outside Dimension
20.12.05	246 x 246 x 50 mm
20.12.07	246 x 246 x 70 mm
20.12.10	246 x 246 x 100 mm



Wire Basket With Lid

Code	Outside Dimension
20.12.05L	246 x 246 x 55 mm
20.12.07L	246 x 246 x 75 mm
20.12.10L	246 x 246 x 105 mm

Middle Size

WIRE BASKETS

Sizes are also compatible with
European-made washer- decontaminators.



Wire Basket

Code	Outside Dimension
20.34.05	405 x 245 x 50 mm
20.34.07	405 x 245 x 70 mm
20.34.10	405 x 245 x 100 mm



Wire Basket With Lid

Code	Outside Dimension
20.34.05L	405 x 245 x 55 mm
20.34.07L	405 x 245 x 75 mm
20.34.10L	405 x 245 x 105 mm

Full Size

WIRE BASKETS

All products comply with DIN58952 range include fine wire mesh baskets for handling small parts.



Wire Basket

Code	Outside Dimension
20.11.05	540 x 245 x 50 mm
20.11.07	540 x 245 x 70 mm
20.11.10	540 x 245 x 100 mm



Wire Basket With Lid

Code	Outside Dimension
20.11.05L	540 x 245 x 55 mm
20.11.07L	540 x 245 x 75 mm
20.11.10L	540 x 245 x 105 mm

Perforated Trays

Wire Baskets and DIN standard baskets with perforated side walls offering greater instrument protection. Satisfies all requirements of instrument sterilization our comprehensive range of stainless steel trays come in a variety of sizes and styles, also available with or without lids.



SHARPLINE.

Half Size

PERFORATED TRAYS



Box With Lid

**Perforated Tray
With Feet**

Code	Outside Dimension
21.12.06L	250 x 259 x 64 mm
21.12.08L	250 x 259 x 84 mm
21.12.11L	250 x 259 x 114 mm



Box


**Perforated Tray
With Feet**

Code	Outside Dimension
21.12.06	250 x 259 x 60 mm
21.12.08	250 x 259 x 80 mm
21.12.11	250 x 259 x 110 mm

Middle Size

PERFORATED TRAYS

A silicone pad is also available to ensure that the wrapped instrument baskets are safely stored on the mesh trays.



Box With Lid

Perforated Tray With Feet

Code	Outside Dimension
21.34.06L	411 x 259 x 64 mm
21.34.08L	411 x 259 x 84 mm
21.34.11L	411 x 259 x 114 mm



Box

Perforated Tray With Feet

Code	Outside Dimension
21.34.06	411 x 259 x 60 mm
21.34.08	411 x 259 x 80 mm
21.34.11	411 x 259 x 110 mm

Full Size

PERFORATED TRAYS

It prevents the baskets from shifting when the mesh rack is transported or moved.



Box With Lid

Perforated Tray With Feet

Code	Outside Dimension
21.11.06L	546 x 259 x 64 mm
21.11.08L	546 x 259 x 84 mm
21.11.11L	546 x 259 x 114 mm



Box


Perforated Tray With Feet

Code	Outside Dimension
21.11.06	546 x 259 x 60 mm
21.11.08	546 x 259 x 80 mm
21.11.11	546 x 259 x 110 mm

Extra Long Size

PERFORATED TRAYS

The composition of the material is checked
by an external knowledge centre.



Box With Lid

**Perforated Tray
With Feet**

Code	Outside Dimension
21.35.06L	671 x 259 x 64 mm
21.35.08L	671 x 259 x 84 mm
21.35.11L	671 x 259 x 114 mm



Box

**Perforated Tray
With Feet**

Code	Outside Dimension
21.35.06	671 x 259 x 60 mm
21.35.08	671 x 259 x 80 mm
21.35.11	671 x 259 x 110 mm

Extra Large Size

PERFORATED TRAYS

The baskets are compatible with all types of detergents, including strong alkalis.



Box With Lid

**Perforated Tray
With Feet**

Code	Outside Dimension
21.36.06L	539 x 346 x 64 mm
21.36.08L	539 x 346 x 84 mm
21.36.11L	539 x 346 x 114 mm



Box

**Perforated Tray
With Feet**

Code	Outside Dimension
21.36.06	539 x 346 x 60 mm
21.36.08	539 x 346 x 80 mm
21.36.11	539 x 346 x 110 mm

STERILIZATION CONTAINERS

Accessories

Comprehensive SHARPLINE accessories for multiple applications and to protect your sensitive instruments.



SHARPLINE.

Accessories


LABELS - SILICONE HANDLE -
IDENTIFICATION LABELS

LABELS



50.01

100 pcs. With indicator.
For Small size and
special containers





50.02.

100 pcs. With indicator.
For Full size and
Half size containers





SILICONE HANDLE & IDENTIFICATION LABELS



Coloured identification labels engraved to your
specifications maximum 16 letters.

Silicone handle has available in six different
colors and detachable by user.



Accessories

SECURITY SEALS & FILTERS

SECURITY SEAL



Code 52.01

100 pcs. with indicator.
made of plastic



Code 52.02

100 pcs. made of plastic



FILTERS (for Standard & Plasma model)

Paper Filter



Code 53.01

Diameter 184 mm

Single use 100 pcs
with indicator.

PTFE(Teflon) Filter




Code 53.02

Diameter 194 mm

Can be applied 5600
cycle times.

Polypropylene Filter



Code 53.03

Diameter 184 mm

Single use 100 pcs.
Diameter: 184 mm
**for hydrogen peroxide
plasma sterilization**



Accessories

FILTERS

FILTERS (For Quarter size, 1/2 Quarter Size & Mini Size)

Paper Filter



Code 53.04

Diameter 165 X 235 mm

Single use 100 pcs
with indicator.

PTFE(Teflon) Filter



Code 53.05

Diameter 150 X 223 mm

Number of sterilization
cycle > 5600



Code 53.06

Diameter 130 X 130 mm

Single use 100 pcs
with indicator.



Code 53.07

Diameter 130 X 130 mm

Number of sterilization
cycle > 5600



Code 53.08

Diameter 230 X 120 mm

Single use 100 pcs
with indicator.



Code 53.09

Diameter 215 X 100 mm

Number of sterilization
cycle > 5600



Accessories


SILICONE PAD

SILICONE PAD

Code 54.01
270 x 125 mm



Code 54.02
270 x 125 mm




Code 54.03
245 x 245 mm




Code 54.04
245 x 245 mm




Code 54.05
225 x 235 mm




Code 54.06
370 x 235 mm



Code 54.07
450 x 235 mm




Code 54.08
500 x 235 mm



Code 54.11
270 x 125 mm



Code 54.12
270 x 125 mm



Code 54.13
245 x 245 mm




Code 54.14
245 x 245 mm



Code 54.15
225 x 235 mm





Code 54.16
370 x 235 mm



Code 54.17
450 x 235 mm




Code 54.18
500 x 235 mm




Can be cut to the
desired dimensions

Code 54.10
545 x 255 mm




Accessories

RACKS



RACKS



Rack, for endoscopy

Holds 10 endoscopic instruments
 Fits into a 1/1 size 20 cm deep container
 Wire basket 1/1 50 mm
 Container 1/1 200 mm
 Silicone pad 520 x 230 mm

Code	Outside Dimension
40.11.01	505 x 250 x 125 mm



Rack, for endoscopy

Holds 10 endoscopic instruments
 Fits into a Extra Long size 18 cm deep container
 Wire basket Extra Long 50 mm
 Container Extra Long 180 mm
 Silicone mat 340x230 (2pcs)



Code	Outside Dimension
40.35.01	550 x 250 x 130 mm



Accessories

RACKS

RACKS



Rack, for arthroscopy

Holds 20 arthroscopy punches
Fits into a 1/1
size 20 cm deep container
Wire basket 1/1 50 mm
Container 1/1 200 mm
Silicone pad 440 x 230 mm

Code	Outside Dimension
40.11.02	515 x 250 x 165 mm

Accessories

RACKS



Rack, for kerrison

Holds 10 kerrison punches
Fits into a 1/2, 1/1 20 cm and Extra
Long size 18 cm deep container
Wire basket 1/2 50 mm
Wire basket 1/1 50 mm
Wire basket Extra long 50 mm
Container 1/2 200 mm
Container 1/1 200 mm
Container Extra long 180 mm

Code	Outside Dimension
40.12.03	220 x 250 x 160 mm

SHARPLINE.

[sharpline.de](https://www.sharpline.de)